

## Printing Labsafe barcodes on TSC TX300 printer

### The directory where the required files is:

|                 |
|-----------------|
| C:\TEMP\PRINTER |
|-----------------|

There is a batch file, called *run.bat* that has to be run for the printer to process jobs. Once the file running – the command line window will appear. The information about the barcodes being processed will be printed in that window.

### Troubleshooting:

1. When running *run.bat*, the command window appears and disappears – this means that the configuration in the *run.bat* files is invalid. Review the batch file content.
2. *run.bat* is running and printer is on but nothing is printed.
  - a. Check your default printing queue in ‘Settings’ under ‘Logged in as: ...’ in the right top corner of Labsafe. The default printing queue name should match the one that is in the batch file (the one that the printer listens to).
  - b. If point *a* is not applicable - restart the printer.
3. Printer produces a number of blanc labels. This is caused by the incorrect printer settings. Run the *DiagTool* and check that the paper width and height correspond to the actual label size.
  - a. If you have previously saved settings for the printer – *Load* them and *Set* using the corresponding buttons in the *DiagTool*. You may need to press the ‘Feed’ button on the printer to set the labels into the right position.
  - b. If you do not have pre-saved settings. Set the paper width manually and use the ‘Calibrate sensor’ button and then use ‘Calibrate’ button under the ‘Auto calibration’ section. The printer then will feed through 3 labels and update its settings. You can use the ‘Get’ button to read the updated settings and then use ‘Save’ button to store them in a file.
4. The labels are not vertically aligned. Calibrate the sensor as described in step *3b*.